

## Product datasheet for **RR204709L3V**

### Men1 (NM\_019208) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Men1 (NM_019208) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Men1
Synonyms:	menin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019208
ORF Size:	1830 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204709).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019208.3</a> , <a href="#">NP_062081.1</a>
RefSeq Size:	2996 bp
RefSeq ORF:	1833 bp
Locus ID:	29417
UniProt ID:	<a href="#">Q9WVR8</a>
Cytogenetics:	1q43
Gene Summary:	displays conservation with the human and mouse homologs for amino acid residues that are mutated in human multiple endocrine neoplasia type 1 [RGD, Feb 2006]



[View online »](#)