

## Product datasheet for RR204717L4

### Pnp (NM\_001106031) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pnp (NM_001106031) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Pnp
Synonyms:	Np
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204717).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



**Locus ID:** 290029

**Cytogenetics:** 15p14

**Gene Summary:** mouse purine-nucleoside phosphorylase (Pnp) exhibits purine-nucleoside phosphorylase and glycosyl transferase activities [RGD, Feb 2006]