

Product datasheet for **RR204722L4V**

Pir (NM_001009474) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pir (NM_001009474) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pir
Synonyms:	MGC109484
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009474
ORF Size:	873 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204722).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009474.1 , NP_001009474.1
RefSeq Size:	1537 bp
RefSeq ORF:	876 bp
Locus ID:	363465
UniProt ID:	Q5M827
Cytogenetics:	Xq14



[View online »](#)

Gene Summary:

Transcriptional coregulator of NF-kappa-B which facilitates binding of NF-kappa-B proteins to target kappa-B genes in a redox-state-dependent manner. May be required for efficient terminal myeloid maturation of hematopoietic cells. Has quercetin 2,3-dioxygenase activity (in vitro) (By similarity).[UniProtKB/Swiss-Prot Function]