

Product datasheet for **RR204743L4V**

Rgs14 (NM_053764) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rgs14 (NM_053764) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rgs14
Synonyms:	MGC108631
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053764
ORF Size:	1632 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204743).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053764.1 , NP_446216.1
RefSeq Size:	2854 bp
RefSeq ORF:	1635 bp
Locus ID:	114705
UniProt ID:	O08773
Cytogenetics:	17p14
Gene Summary:	binds Galpha(i1)GDP and acts as a GTPase activator by maintaining G-protein subunit dissociation; regulates G-protein mediated signaling [RGD, Feb 2006]



[View online »](#)