

## Product datasheet for **RR204804L3V**

### Apoc2 (NM\_001085352) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Apoc2 (NM_001085352) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Apoc2
Synonyms:	RGD1560725
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001085352
ORF Size:	291 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204804).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001085352.1</a> , <a href="#">NP_001078821.1</a>
RefSeq Size:	518 bp
RefSeq ORF:	294 bp
Locus ID:	292697
Cytogenetics:	1q21
Gene Summary:	cofactor and activator of lipoprotein lipase; crucial for the hydrolysis of triacylglycerols and very-low-density lipoproteins; mRNA found mainly in the liver [RGD, Feb 2006]



[View online »](#)