

Product datasheet for **RR204827L3V**

Pdhb (NM_001007620) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdhb (NM_001007620) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pdhb
Synonyms:	MGC94149
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007620
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204827).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007620.1 , NP_001007621.1
RefSeq Size:	1816 bp
RefSeq ORF:	1080 bp
Locus ID:	289950
UniProt ID:	P49432
Cytogenetics:	15p14
Gene Summary:	The pyruvate dehydrogenase complex catalyzes the overall conversion of pyruvate to acetyl-CoA and CO ₂ , and thereby links the glycolytic pathway to the tricarboxylic cycle. [UniProtKB/Swiss-Prot Function]



[View online »](#)