

Product datasheet for **RR204830L4V**

Tgfa (NM_012671) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tgfa (NM_012671) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tgfa
Synonyms:	RATTGFAA; TGFAA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012671
ORF Size:	477 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204830).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012671.2 , NP_036803.1
RefSeq Size:	4222 bp
RefSeq ORF:	480 bp
Locus ID:	24827
UniProt ID:	P01134
Cytogenetics:	4q34
Gene Summary:	acts as an EGF receptor ligand; inhibits apoptosis that results from epithelial detachment from the extracellular matrix [RGD, Feb 2006]



[View online »](#)