

Product datasheet for **RR204836L3V**

Acad11 (NM_001108181) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acad11 (NM_001108181) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Acad11
Synonyms:	RGD1306270
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001108181
ORF Size:	2337 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204836).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001108181.1 , NP_001101651.1
RefSeq Size:	3449 bp
RefSeq ORF:	2340 bp
Locus ID:	315973
UniProt ID:	B3DMA2
Cytogenetics:	8q32
Gene Summary:	Acyl-CoA dehydrogenase, that exhibits maximal activity towards saturated C22-CoA. [UniProtKB/Swiss-Prot Function]



[View online »](#)