

## Product datasheet for RR204843L4V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Tyro3 (NM\_017092) Rat Tagged ORF Clone Lentiviral Particle

## **Product data:**

**Product Type:** Lentiviral Particles

**Product Name:** Tyro3 (NM\_017092) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Tyro3
Synonyms: Brt

Mammalian Cell Puromycin

Selection:

. a. o. . . y c. . .

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

**ACCN:** NM\_017092 **ORF Size:** 2640 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RR204843).

Sequence:
OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**RefSeq:** <u>NM 017092.1, NP 058788.1</u>

 RefSeq Size:
 3726 bp

 RefSeq ORF:
 2643 bp

 Locus ID:
 25232

 UniProt ID:
 P55146

 Cytogenetics:
 3q35

**Gene Summary:** receptor tyrosine kinase; may have a role in development or function of specific neuronal

populations [RGD, Feb 2006]

