

## Product datasheet for **RR204886L4V**

### Hbp1 (NM\_013221) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hbp1 (NM_013221) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hbp1
Synonyms:	Hmgb1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013221
ORF Size:	1569 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204886).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013221.2</a> , <a href="#">NP_037353.2</a>
RefSeq Size:	2671 bp
RefSeq ORF:	1572 bp
Locus ID:	27080
UniProt ID:	<a href="#">Q62661</a>
Cytogenetics:	6q16
Gene Summary:	a transcriptional regulator that may play a central role in the stabilization and/or assembly of several multifunctional complexes through protein-protein interactions [RGD, Feb 2006]



[View online »](#)