

## Product datasheet for **RR204892L4V**

### **Bod1 (NM\_001013854) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bod1 (NM_001013854) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bod1
Synonyms:	Fam44b; RGD1305622
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001013854
ORF Size:	711 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204892).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001013854.1</a> , <a href="#">NP_001013876.2</a>
RefSeq Size:	1486 bp
RefSeq ORF:	714 bp
Locus ID:	287173
UniProt ID:	<a href="#">Q6AYJ2</a>
Cytogenetics:	10q12
Gene Summary:	Required for proper chromosome biorientation through the detection or correction of syntelic attachments in mitotic spindles.[UniProtKB/Swiss-Prot Function]



[View online »](#)