

## Product datasheet for **RR204893L4V**

### Thop1 (NM\_172075) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Thop1 (NM_172075) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Thop1
Synonyms:	EP24.15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172075
ORF Size:	2061 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204893).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172075.2</a> , <a href="#">NP_742072.2</a>
RefSeq Size:	2408 bp
RefSeq ORF:	2064 bp
Locus ID:	64517
UniProt ID:	<a href="#">P24155</a>
Cytogenetics:	7q11
Gene Summary:	metallopeptidase that catalyzes the thiol-dependent cleavage of many neuropeptides [RGD, Feb 2006]



[View online »](#)