

Product datasheet for **RR204900L3V**

Ctrc (NM_001077649) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctrc (NM_001077649) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ctrc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001077649
ORF Size:	804 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204900).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001077649.2 , NP_001071117.1
RefSeq Size:	898 bp
RefSeq ORF:	807 bp
Locus ID:	362653
UniProt ID:	P55091
Cytogenetics:	5q36
Gene Summary:	Has chymotrypsin-type protease activity and hypocalcemic activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)