

Product datasheet for **RR204903L4V**

Eaf2 (NM_172047) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eaf2 (NM_172047) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Eaf2
Synonyms:	Traits
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172047
ORF Size:	786 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204903).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172047.1 , NP_742044.1
RefSeq Size:	1132 bp
RefSeq ORF:	789 bp
Locus ID:	266787
UniProt ID:	Q811X5
Cytogenetics:	11q21
Gene Summary:	Acts as a transcriptional transactivator of ELL, ELL2 and TCEA1 elongation activities (By similarity). Potent inducer of apoptosis in prostatic and non-prostatic cell lines. [UniProtKB/Swiss-Prot Function]



[View online »](#)