

Product datasheet for **RR204944**

Wwp1 (NM_001024757) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wwp1 (NM_001024757) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wwp1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide
Sequence:**

>RR204944 representing NM_001024757
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCACTGCTTACCAAGATCTGACACTAGTGATAACCCACAGTGAAGGTTACAGTTAAAGGTAACAG
 TTTCTAGTGCCAACTAAAAAGAAAAAGAACTGGTTTGAACAGCAATATATACTGAAGTCATAGTAGA
 TGGAGAAATTAAGAAAACAGCAAAATCCAGTAGTCTTCTAATCCAAAATGGGATGAGCAGCTGATTGTA
 AATGTGGCCACAGACCACCTTGGAGTTCGAGTTTGGAGTCATCACACTTTAAAGGCAGACGCTCTAT
 TAGGAAAAGCAACAGTGGACTTGAACAAGTGCTGTTAACACACAACAGAAAATTAGAAAAAGTAAAGA
 ACAACTAAAACCTTCTTGGAAAACAAGAAATGGAATAGTACAACTGGGAACTGACAGTTGCTCTTGAT
 GGATTAGTGATTGAGCCAGAAGCTGTGACAAACCGCAGCTCCTCGCCACCCATAGAAATTCAGCAGAATG
 GGGATGCTGTCCATGAGAATGGAGACCCTTTGACAAGGACAACCTCCAGGTTGGCTGTTGAAGGACTAT
 TGGAAATGATAATCATGTATCAACAAGTACTGTTGTTCCGAGCTCATGCTGTTACATATAGTTAATGGA
 GAAAACACACCTTCATCCCCGCCTCAGGTTGCTACCAGACCCAAAATACTCCAGCTCCAAAACCACTCA
 CGTCTGAGCCACCAGTGACACAGTTAATGGAGAATCATCATCTGTCTAGCTGATAGTACTTCAGCCAT
 GGGTACTTCACTACCATCTGAAGACTCCACCTCGACTTCCAATTGCACTAGCACCACCGTCCAAGAGCCT
 CCTGTACAAGAACCCCGCAGTCCCTCAGAACACAGTGAATGTAATCCCTCTGCCAGTGCAGAAGTGGGAC
 CAGAAGCTCGGAGCCTCATAGACCCTGACTCTGATTCTAGAAAATAATTCTGGTTTTGACAAAAGTAAAGCA
 GCCAGAAGGGTGTGTGGAACCTCTCCGGCCACAGTCGGGAAATACCAACACAGAATCTTTGCCATCAGGG
 TGGGAACAGAGAAAAGATCCTCATGGTAGAACCTATTATGTGGATCATAACTCGAACTACTACATGGG
 AAAGACCACAACCTTTGCCCCAGGTTGGGAAAGAAGGTTGATGATCGTGAAGAGTTTATTATGTGGA
 TCATAACACCAGAACAACAACCTGGCAGCGACCTACCATGGAGTCTGTTCCGGAATTTGAACAGTGGCAG
 TCTCAGCGGAATCAATTGCAGGGAGCTATGCAACAATCAACCAACGATACCTCTATTCGGCTTCAATGT
 TAGCTGCAGAAAATGACCCATATGGACCCTTACCACCAGGCTGGGAAAAAGAGTGGATTCAACAGACAG
 AGTTTACTTTGTGAATCATAATACAAAAACAACCCAGTGGGAAGATCCAAGAAGTCAAGGTTTACCAAT
 GAAGAACCCTCCAGAAGGCTGGGAAATTAGGTATACTCGGAAAGGTGTCAGGTTTGTGACCATA
 ACACAAGAACAACAACATTCAAAGATCCCCGCAATGGGAAGTATCGGTAACAAAGGTGGTCCGCAGAT
 CGTTATGAGCGCAGCTTAGGTGGAACTGGCCACTTCGTTACTTGTGCCAGTCAATGCACTACCC
 AGTCATGTAAGATCAATGTGTCCGACAGACATTATTTGAAGATTCCTTCAACAGATTATGGCATTAA
 AACCGTATGACTTGAGGAGGCGATTATATGTAATATTTAGAGGAGAAGAAGGCCTCGATTATGGTGGCCT
 TGCAAGAGAATGGTTTTCTTGCTTTCCATGAAGTTTTGAACCAATGTATTGCTTATTTGAGTATGCA
 GGCAAGAATAACTATTGTTTACAGATAAACCCAGCGTCAACTATAAATCCTGATCATCTTTCATACTTCT
 GTTTTATTGGTCGCTTATTGCAATGGCGCTATCCACGGAAAGTTTATTGATACTGGTTTTCTTTACC
 ATTTCTACAAGCGTATGCTAAGTAAAAAAGTACTATCAAGGATTTGGAATCTATCGATACTGAATTTTAC
 AATTCTCTAATCTGGATAAGAGATAACAATATTGAAGAATGTGGCTTAGAAATGACTTTTCTGTAGACA
 TGGAGATTCTGGGAAAAGTTACTTCACATGATTTAAAGTTGGGAGGTTCCAATATCTGGTACTGAAGA
 AAACAAGATGAATATATTGGTTAATGACAGAATGGCGTTTTTTCACGAGGAGTACAAGAGCAGACCCAAA
 GCTTTCCTTGATGGGTTTAAATGAAGTTGCCCCCTTCAGTGGCTACAGTACTTTGATGAAAAAGAATTGG
 AGGTTATGCTGTGTGGCATGCAAGAGGTTGACCTGGCAGACTGGCAGAGAAAACACTGTTTACCGACATTA
 CACGAGAAAACAGCAAGCAGATTATTTGGTTTTGGCAGTTTGTGAAAGAGACTGATAATGAAGTACGAATG
 CGACTCTGCAGTTCGTCACTGGAACCTGTCGCTTGCCTCTAGGAGGATTTGCTGAGCTCATGGGAAGTA
 ATGGGCCCTCAAAAATTCTGTATTGAAAAAGTTGGCAAAGACACCTGGTTACCAAGAAGCCATACATGTTT
 TAATCGCTTGGATCTACCACCATATAAGAGTTATGAACAACAAAGGAAAAAAGTCTTTTTGCAATAGAA
 GAGACAGAGGGATTTGGACAAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAAGTTTAA

Protein Sequence: >RR204944 representing NM_001024757
 Red=Cloning site Green=Tags(s)

MATASPRSDTSDNHSGRLQLKVTVSSAKLKRKKNWFGTAIYTEVIVDGEIKKTAKSSSSSNPKWDEQLIV
 NVAPQTTLEFRVWSSHHTLKADALLGKATVDLKQVLLTHNRKLEKVKQLKLSLENKNGIVQTGELTVVLD
 GLVIEPEAVTNRSSPPIETIQNGDAVHENGDP LTRTTPLAVEGTIGIDNHVSTSTVVPSSCCSHIVNG
 ENTPSSPPQVATRPQNTAPKPLTSEPTSDTVNGESSVLADSTSAMGTS LPS E D S T S T S N C T S T T V Q E P
 PVQEPPASSEHSECI PSASAEVGP EARS LIDPDS DSRNNSGFDKVRQPEGCV EPLRPQSGNTNTESLPSG
 WEQRKDPHGRTYYVDHNTRTTTWERPQLPPGWERRVDDRGRVYYVDHNTRTTTWQRP T M E S V R N F E Q W Q
 SQRNQLQGAMQQFNQRYLYSASMLAAENDPYG L P P G W E K R V D S T D R V Y F V N H N T K T T Q W E D P R T Q G L P N
 EEPLPEGWEIRYTR E G V R Y F V D H N T R T T T F K D P R N G K S S V T K G G P Q I A Y E R S F R W K L A H F R Y L C Q S N A L P
 SHVKINVS R Q T L F E D S F Q Q I M A L K P Y D L R R R L Y V I F R G E E G L D Y G G L A R E W F F L L S H E V L N P M Y C L F E Y A
 GKNNYCLQINPASTINPDHLSYFCIFGRFIAMALFHGKFIDTGFSLPFYKRMLSKLTIKDESIDTEFY
 NSLIWIRDNNIEECGLEMYFSVDMEILGKVTSHDLKLGGSNILVTEENKDEYIGLMT E W R F S R G V Q E Q T K
 AFLDGFNEVVPLQWLQYFDEKELEVMLCGMQEVDLADWQRNTVYRHYTRNSKQIIWFWQFVKETDNEVRM
 RLLQFVTGTCRLPLGGFAELMGSNGPQKFCIEKVGKDTWLP RSHTCFNRDLDPYKSYEQLKEKLLFAIE
 ETEGFGQE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001024757

ORF Size: 2754 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001024757.1](#), [NP_001019928.1](#)

RefSeq Size: 3331 bp

RefSeq ORF: 2757 bp

Locus ID: 297930

Cytogenetics: 5q13

MW: 104.6 kDa