

## Product datasheet for **RR204964L3V**

### Rpl35 (NM\_212511) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rpl35 (NM_212511) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rpl35
Synonyms:	MGC72958
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_212511
ORF Size:	369 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204964).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_212511.2</a> , <a href="#">NP_997676.1</a>
RefSeq Size:	476 bp
RefSeq ORF:	372 bp
Locus ID:	296709
UniProt ID:	<a href="#">P17078</a>
Cytogenetics:	3q12
Gene Summary:	mouse homolog is Rpl35, a component of the 60S ribosomal protein complex; may be a structural component of the ribosome; may be involved in protein biosynthesis [RGD, Feb 2006]



[View online »](#)