

## Product datasheet for **RR204978L4V**

### Cirbp (NM\_031147) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cirbp (NM_031147) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cirbp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031147
ORF Size:	516 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204978).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031147.2</a> , <a href="#">NP_112409.2</a>
RefSeq Size:	1482 bp
RefSeq ORF:	519 bp
Locus ID:	81825
Cytogenetics:	7q11
Gene Summary:	expression decreased by ischemia and H2O2; may mediate generation of reactive oxygen species [RGD, Feb 2006]



[View online »](#)