

## Product datasheet for **RR204979L3V**

### Tbccd1 (NM\_001012016) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tbccd1 (NM_001012016) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tbccd1
Synonyms:	Crygs
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012016
ORF Size:	1707 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204979).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001012016.1</a> , <a href="#">NP_001012016.1</a>
RefSeq Size:	2167 bp
RefSeq ORF:	1710 bp
Locus ID:	303830
UniProt ID:	<a href="#">Q5FVR8</a>
Cytogenetics:	11q23
Gene Summary:	Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.[UniProtKB/Swiss-Prot Function]



[View online »](#)