

Product datasheet for **RR204988L2V**

Anxa1 (NM_012904) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Anxa1 (NM_012904) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Anxa1
Synonyms:	Anx1; p35
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_012904
ORF Size:	1038 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204988).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012904.2 , NP_037036.1
RefSeq Size:	1401 bp
RefSeq ORF:	1041 bp
Locus ID:	25380
UniProt ID:	P07150
Cytogenetics:	1q51
Gene Summary:	plays a role in regulation of insulin secretion; may inhibit phopsholipase A2 [RGD, Feb 2006]



[View online »](#)