

Product datasheet for RR204988L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Anxa1 (NM_012904) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Anxa1 (NM 012904) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Anxa1

Synonyms: Anx1; p35

Mammalian Cell Puromycin

Selection:

Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_012904 **ORF Size:** 1038 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RR204988).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 012904.2, NP 037036.1

RefSeq Size: 1401 bp
RefSeq ORF: 1041 bp
Locus ID: 25380
UniProt ID: P07150
Cytogenetics: 1q51

Gene Summary: plays a role in regulation of insulin secretion; may inhibit phopsholipase A2 [RGD, Feb 2006]

