

Product datasheet for **RR204990L3V**

Thy1 (NM_012673) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Thy1 (NM_012673) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Thy1
Synonyms:	CD7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012673
ORF Size:	483 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR204990).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012673.2 , NP_036805.1
RefSeq Size:	1328 bp
RefSeq ORF:	486 bp
Locus ID:	24832
UniProt ID:	P01830
Cytogenetics:	8q22
Gene Summary:	cell surface glycoprotein antigen [RGD, Feb 2006]



[View online »](#)