

Product datasheet for **RR205022**

Hoxd10 (NM_001107094) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hoxd10 (NM_001107094) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Hoxd10
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205022 representing NM_001107094
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCCTTTCCCAACAGCTCTCCTGCTGCTAATACTTTTTAGTAGATTCCTTGATCAGTGCCTGCAGGA
GTGACAGTTTTTATTCGAGCAGCGCCAGCATGTACATGCCACCCTAGCGCAGACATGGGACCTATGG
AATGCAAACCTGTGGACTGCTCCCGTCTCTGGCCAAAAGAGAAGTGAACCACAAAATATGGGTATGAAT
GTACATCCTTATATACCTCAAGTAGACAGTTGGACAGACCCGAACAGATCTTGCGAATAGAGCAACCTG
TTACACAGCAAGTCCCCACTTGCTCATTACCCGCAACATTAAGAAGAATCCAATTGCTGCATGTATTC
TGATAAGCGCAACAACTCATTCTGCCGAGGTCCTTGTACCAGAGGCTGGTCCCGAGTCTTGTCTCT
GTTGAGAACCCCGAGGTTCTGTCCCTGGATATTTAGACTGAGTCAAGCCTACGCCACCGGAAAACCC
AAGAGTACAATAACAACCCGAAGGCAGCTCCACTGTATGCTCCAGCTCAACCCCTGTGGCGGGCCAA
GCCGAGTCTCGGCTGCCAGCTGCAGATGGAAAAGAAGATGAACGAAACCGCGAGCGGCCAGGAGCCC
ACTAAAGTCTCCAAGTGGAGAGCCCCGAGGCCAAAGGCGGCCTCCAGAAGAGAGGAGCTGCCTGGCTG
AGGTCTCCGTGTCCAGCCCCGAAGTGCAGGAGAAGGAAAGCAAAGAGGAAATCAAGTCTGATACTCCAAC
CAGCAATTGGCTCACTGCAAAGAGTGGCAGAAAAGAAGAGGTGCCCTTACACTAAACACCAACGCTGGAA
TTAGAAAAAGAGTTCTTGTTCATATGTACCTCACCCGCGAGCGCCGCTAGAGATCAGTAAAGAGCGTTA
ACCTCACCGACAGGCAGGTCAAGATTTGGTTTCAAACCAGCAATGAACTCAAGAAGATGAGCCGAGA
GAATCGAATCCGAGAAGTACCGCAACCTCACCTTTTCT

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR205022 representing NM_001107094
 Red=Cloning site Green=Tags(s)

MSFPNSSPAANTFLVDSLISACRSDSFYSSASMYPPPSADMGTYGMQTCGLLPSLAKREVNHQNMGMN
 VHPYIPQVDSWTDPNRSCRIEQVPTQQVPTCSFTANIKESNCCMYSDKRNLISAEVPSYQRLVPESCP
 VENPEVPVPGYFRLSQTATGKTQEYNNNPEGSSTVMLQLNPRGAAPQLSAAQLQMEKKMNETASGQEP
 TKVSQVESPEAKGGLPEERSCLAEVSVSSPEVQEKEKEEIKSDTPTSNWLTAKSGRKKRCPYTKHQTL
 LEKEFLFNMYLTRERRLEISKSVNLTDRQVKIWFQNRMMKLLKMSRENRIRELTANLTFSS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107094

ORF Size:	1020 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107094.1 , NP_001100564.1
RefSeq Size:	1848 bp
RefSeq ORF:	1023 bp
Locus ID:	303991
Cytogenetics:	3q23
MW:	38.4 kDa