

## Product datasheet for **RR205035L3V**

### Fcer1g (NM\_001131001) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Fcer1g (NM_001131001) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fcer1g
Synonyms:	MGC188226
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001131001
ORF Size:	258 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205035).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001131001.1</a> , <a href="#">NP_001124473.1</a>
RefSeq Size:	657 bp
RefSeq ORF:	261 bp
Locus ID:	25441
UniProt ID:	<a href="#">P20411</a>
Cytogenetics:	13q24
Gene Summary:	subunit of high-affinity receptor for immunoglobulin E; mediates the release of mediators responsible for allergic symptoms [RGD, Feb 2006]



[View online »](#)