

Product datasheet for **RR205067**

F11 (NM_001047848) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	F11 (NM_001047848) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	F11
Synonyms:	LRRGT00086
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR205067 representing NM_001047848
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGATCCCCAGGATTAACCTTACATCATTACTCTTGGTGCCTCGCCAGCGTGAAACTTACTCACAA
 GAACCAAGTTAGGGAGTCTGCCCATGTGCCATCAGTAGACAAATGCATAAAGCAAATGCTGGAAGTACG
 CACAAAGTTCCGCCCTTGTACTTCTGTCCCGCTATTGGCTGTCAGATCTTTACTATAACCAATCAGCA
 GTAAGACAACGCTGGAACATCTCTTTGGTTGCCTGTTGGCAATATGGGCTGTTGGGAGAACAGGTTTG
 CCTGCATCCTGAAGGACAGTCTCACAGAAACATTGCCGATGGTACACATGACAGGAGCGGTTTCTGGATA
 TTCCTTCAAGCAATGCCCTCAGCAATTAAGTGCTTGCAGCAAAAACGTGTACGTGAACCTGGACATGAAG
 GGCATGAACATAACAGTTCTGCCGTGAAGAATGCTCGAGAATGTCAGGAGAGATGCACCAACGATGCC
 ACTGCCAGTTTTTACATACGCAACGGGGCATTTCAGTGTGGACCATCGTTTCTCTCATCGTACAGC
 TTGTATCAGGGACATTTCCCGAACCCGTGCTTGCAGACCTAACATTGACAGCGTGTGGCCCCAGAT
 GCTTTTGTCTGTCTGCGATTTGCACGCATCACCCACGTGCCTGTTTTTACATTCTTTACCCAAGCGT
 GGCCAAAAGAATCCCAGAGACATCTTGGCTCCTTAAAACATCTGAAAGCGGATTTCCAAGCACACGCAT
 CACAAAGAGCAACGCCCTTCCGGGTTCACTCCAGCACTGCAGGCACAGCATCCCAAGTATTCTGTCTAT
 CCGAACTTTTACAACGACTGATTTCTTGGGAGAAGAACTGGACATCGTCGATGTGAAAGGCAAAGAAT
 CCTGCCAGAAAATGTGCTCTGACAACGTCGCTGCCAGTCTTTACCTATTATCCCTCTCGAGGATCTTG
 CAGTGAGAGGAATATGCACAACCAAAATCAGGCCACGGGTGTTCCGAGGAGCTGCGTCTGTTACGGCGA
 GTGGCCATGGCAGGTGACCTGCACACCACCCAGGGACACCTGTGTGGAGGCTCCATCATTGAAACCGG
 TGGATATTGACAGCGGCTCATTGTTTCTCTGGTTAGTATTAACGATATTCGCCTGGGTCGTCAACC
 AAGTACAGAGACACCTAAAACCTCTGCGTGTCTACGGTGGTATTGTAATCAATCAGAAATAAATGAAGAT
 ACCACTTTCTTCAAGGTTCAAGAAATGATAATTATGATCAATATACATCGGCAGAAAAGTGGTTTGACA
 TTGCCCTCTTAAAACCTGGAACCGGCCATGAATTACACAGATGAAGTACAAAGTACTCTCCAGAAAGCCAA
 GGTACCATTGGTGTGCAATGAAGAATGTCAAACAAGATACAGAAAACATAAAAATAACCAACAAGGTGATC
 TGTGCAGGATATAAGGAAGGAGGAAAGGATACGTGTAAGGGAGATTCTGGAGGGCCCTGTCTGCAAAAC
 ACAATGGGGTCTGGCACTTGGTGGGCATACAAGCTGGGGTGAAGGCTGCGGCCAGAAAGAGAGGCCGGG
 TGCTACACCAACGTGGCCAAGTATGTGGACTGGATTTGGAGAAAACCTAGTCGGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR205067 representing NM_001047848
 Red=Cloning site Green=Tags(s)


MEIPRIKLTSLLLVPRQRGNLLTRTKLGS LPMCP SVDKCIKQMLEVTRTKFRPCTSCPGYWLSDLYNQSA
 VRQRWNISLVACWQYGLFGQNR FACILKDSL TETLPMVHMTGAVSGYSFKQCPQQLSACSKNVYVNLDMK
 GMNYNSSAVKNARECQERCTNDAH CQFFTYATGHFSPVDHRF SHRTACIRDIFPNTVLADLNIDSVLAPD
 AFVCRRICTHPTCLFFFTTQAWPKESQRHLCLLKTSESGFPSTRITKSNALSGFSLQCHRHSIPVFCH
 PNFYNDTDFLGEELDIVDVKGKESCQKMCSDNVR CQFFTYPSRGSCSERNMHNQNAQQVRRSCVCSRR
 VAMAGDPAHHPGTPVWRLHWHKPV DIDSGLFLWVSINDIRLGS SSPSTETPKTLRVYGGIVNQSEINED
 TTFFRVQEMI IHDQY TSAESGFDI ALLKLEPAMNYTDEVQSTLQKAKVPLVSNEECQTRYRKHK I TNKVI
 CAGYKEGGKDTCKGDSGGPLSCKHNGVWHL VGITSWGEGCGQKERPGVYTNVAKYVDWILEKTQSE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001047848
ORF Size:	1668 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001047848.1 , NP_001041313.1
RefSeq Size:	1671 bp
RefSeq ORF:	1671 bp
Locus ID:	290757
Cytogenetics:	16q11
MW:	63 kDa