

Product datasheet for **RR205097**

Eif3k (NM_001106242) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eif3k (NM_001106242) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Eif3k
Synonyms: Eif3s12; PLAC-24
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR205097 representing NM_001106242
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGAGCCATGGAGGGATCAAGCCCTGGATGAGAACATGGACCTCCTGGAAGGCATAACTGGCTTTG
AAGACTCTGTCCGAAAATTTATCTGCCATGTCGTGGGCATCACCTACCAGCACATCGACCGCTGGCTGCT
GGCCGAGATGCTCGGAGACCTGACTGACAACCAGCTGAGGGTGTGGATGAGCAAGTACGGCTGGAGCGCG
GACGAGTCAGGGCAGGTCTTCATCTGCAGCCAGGAGAGCATTAAAGCCCAAGAACATCGTGGAGAAGA
TCGACTTTGACAGTGTGCCAGCATCATGGCCTCGTCCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR205097 representing NM_001106242
Red=Cloning site **Green**=Tags(s)

MQRPWRDQALDENMDLLEGITGFEDSVRKFICHVVGITYQHIDRWLLAEMLDLTDNQLRVWMSKYGWSA
DESGQVFICSQEESIKPKNIVEKIDFDSVSSIMASSQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-Mlul



[View online »](#)

Cloning Scheme:


ACCN: NM_001106242

ORF Size: 321 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001106242.1](#), [NP_001099712.1](#)

RefSeq Size: 481 bp

RefSeq ORF: 324 bp

Locus ID: 292762
 Cytogenetics: 1q21
 MW: 12.4 kDa

Product images:



Circular map for RR205097