

Product datasheet for **RR205106**

Dph1 (NM_001105809) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dph1 (NM_001105809) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dph1
Synonyms: RGD1562694
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205106 representing NM_001105809
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGGCCTGGTTGTGCCGAGACTCCGAGCCAGGAAGCCGAGTCGGCCCTGGCAGAGGTGCGATCT
 CTCGGGGCCGACTGGCCAATCAGATCCCCCTGAGATCCTGAACAGTCCCAGCTACAGGCCGCTGTCCA
 TGCCCTGCCTTAACTACAACCTTGAGATCCCCAAGACCATCTGGAGGATCCAGCAAGCCAGGCCAAG
 AAGGTGGCCTTACAAATGCCAGAAGGCCTCTGCTCTTTGCTGCACCATTGTGGATATCTTGGAAAGGT
 TCACAAAGGCTGAGGTGATGGTATGGCGATGTGACCTACGGAGCATGCTGTGTGGACGACTTCACTGC
 AAGGGCCTTGGGAGTTGACTTCCCTGGTGCATATGGCCACAGCTGCCTAGTCCCCATGGACACCTCAGTC
 CAAGACTTCCGTGTGCTGTATGTCTTTGTGGATATCCGGATAGACACTGCCACCTTCTGGACTCGGTCC
 GCCTCACCTTACCCAGGCAGCTCGCTTGCTCTGGTCAGCACCATTCAGTTTGTGTCAACTTACAGGC
 AGCTGCCAGGAGCTGAAAGCTGATTATCACATCAGTGTCCACAGTGAAGCCCTGTCCCCTGGGGAG
 ATCCTAGGCTGCACGTCCCCTCGACTACCAAGGAAGTGAAGCTGTTGTGTATCTTGGAGATGGCCGCT
 TCCATCTGGAGTCTGTCATGATCGCAACCCTAATATACCTAGTTACCGGTACGACCCATATAGCAAAGT
 CCTATCCAGAGAATACTATGACCATCAGCGCATGCAGGCCACTCGTCAGGAAGCCATCGCTGCTGCACGC
 TCAGCCAAGTTCTGGGCCTTATTCTGGGAACCTTTGGGCCCGCAGGGAAGTCCCAAGTCTGGAGCACT
 TGGAACTCAGCTCAGAACTTGGGACTTCTTTCTGAGGCTGCTTCTCTGAGATCTTCCCCAGCAA
 GCTCAGTCTACTTCCCTCGGTGGACGTGTGGTGCAGGTGGCATGTCCACGCCTTCCATTGACTGGGGC
 TCAGCCTTTCCCAAGCCACTGCTGACACCCTACGAGGCAGCTGTGGCCCTGAAAGAGATTTCTTGGCAGC
 AACCTACCTATGGACTTCTACGCTGGCAGCTCCTTAGGGCCATGGACAGTGAACCATGGTCGGGACCG
 AGCACCCAGGGTCTCTGCCAGCCTGCATCCGACAAGGTGCAGCAGGGTCCAGAGGCCACTCTCCAGTC
 CCGGCCTGTGAGGGCTGCAGCTGCGCAGACCAGAAAGCTACACCGCCAGCTCCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR205106 representing NM_001105809
 Red=Cloning site Green=Tags(s)


MAALVVPETSEPGSRVGPGRGRISRGLANQIPPEILNSPQLQAAVHALPSNYNFEIPKTIWRIQQAQAK
 KVALQMPEGLLLFACTIVDILERFTKAEVMVMGDVTYGACCVDDFTARALGVDFLVHYGHSLVPMDSV
 QDFRVLVYFVDIRIDTAHLLDSVRLTFTPGSSLALVSTIQFVSTLQAAAQELKADYHISVPQCKPLSPGE
 ILGCTSPRLPKVEAVVYLGDRFHLESVMIANPNIPSYRYDPYSKVL SREYYDHQRMQATRQEIAAAR
 SAKFWGLILGTLGRQSPKVLHLESQRLNLGLPFLRLLSEIFPSKLSLLPSVDVWVQVACPRLSIDWG
 SAFFPKLLTPYEAAVALKEISWQQPYPMDFYAGSSLPWTVNHGRDRAPRGLCQPASDKVQQGSRGHSPV
 PACEGCSADQKATPPAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001105809

ORF Size:	1314 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001105809.2 , NP_001099279.1
RefSeq Size:	2075 bp
RefSeq ORF:	1317 bp
Locus ID:	287523
Cytogenetics:	10q24
MW:	48.1 kDa