

Product datasheet for **RR205126L4V**

Dpysl3 (NM_012934) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpysl3 (NM_012934) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dpysl3
Synonyms:	Crmp4; TUC-4b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012934
ORF Size:	1710 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205126).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012934.1 , NP_037066.1
RefSeq Size:	1792 bp
RefSeq ORF:	1713 bp
Locus ID:	25418
UniProt ID:	Q62952
Cytogenetics:	18p11-q11
Gene Summary:	protein expressed in the nervous system during development that may be involved in neuronal plasticity, neurite growth and extension, and may serve as a ligand of the extracellular matrix during brain development [RGD, Feb 2006]



[View online »](#)