

## Product datasheet for **RR205142L3V**

### Trpv2 (NM\_017207) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Trpv2 (NM_017207) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Trpv2
Synonyms:	Vrl1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017207
ORF Size:	2283 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205142).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017207.3</a> , <a href="#">NP_058903.2</a>
RefSeq Size:	2834 bp
RefSeq ORF:	2286 bp
Locus ID:	29465
UniProt ID:	<a href="#">Q9WUD2</a>
Cytogenetics:	10q23
Gene Summary:	ion channel that may be activated by high temperature; may transduce high threshold heat responses in sensory ganglia [RGD, Feb 2006]



[View online »](#)