

Product datasheet for **RR205144L4V**

Yme111 (NM_053682) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Yme111 (NM_053682) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Yme111
Synonyms:	FtsH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053682
ORF Size:	2145 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205144).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053682.2 , NP_446134.2
RefSeq Size:	3006 bp
RefSeq ORF:	2148 bp
Locus ID:	114217
Cytogenetics:	17q12.3
Gene Summary:	human homolog has similarity to mitochondrial AAA proteases, especially to yeast Yme1p; may be a candidate gene for a form of hereditary spastic paraplegia or neurodegenerative disorder [RGD, Feb 2006]



[View online »](#)