

## Product datasheet for **RR205165L3V**

### Itgal (NM\_001033998) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Itgal (NM_001033998) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Itgal
Synonyms:	LFA-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033998
ORF Size:	3480 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205165).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001033998.2</a> , <a href="#">NP_001029170.1</a>
RefSeq Size:	4363 bp
RefSeq ORF:	3483 bp
Locus ID:	308995
Cytogenetics:	1q37
Gene Summary:	human homolog encodes the integrin alpha L chain which combines with beta 2 chain to form LFA-1 [RGD, Feb 2006]



[View online »](#)