

Product datasheet for **RR205178L3V**

Kdelr1 (NM_001017385) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kdelr1 (NM_001017385) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kdelr1
Synonyms:	MGC109169
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001017385
ORF Size:	636 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205178).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001017385.1 , NP_001017385.1
RefSeq Size:	1476 bp
RefSeq ORF:	639 bp
Locus ID:	361577
UniProt ID:	Q569A6
Cytogenetics:	1q22
Gene Summary:	Receptor for the C-terminal sequence motif K-D-E-L that is present on endoplasmic reticulum resident proteins and that mediates their recycling from the Golgi back to the endoplasmic reticulum.[UniProtKB/Swiss-Prot Function]



[View online »](#)