

Product datasheet for **RR205220L3V**

Map4k3 (NM_133407) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Map4k3 (NM_133407) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Map4k3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133407
ORF Size:	2619 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205220).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133407.1 , NP_596898.1
RefSeq Size:	4109 bp
RefSeq ORF:	2622 bp
Locus ID:	170920
UniProt ID:	Q924I2
Cytogenetics:	6q11
Gene Summary:	a mitogen-activated protein kinase kinase kinase kinase that functions upstream of c-Jun N-terminal kinase [RGD, Feb 2006]



[View online »](#)