

Product datasheet for **RR205234**

Dnm1l (NM_053655) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dnm1l (NM_053655) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dnm1l
Synonyms:	DLP1; Dnml1; Drp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR205234 representing NM_053655
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGCGCTGATCCCGTTCATCAATAAGCTGCAAGACGTCTTCAACACCGTGGGCGCCGACATCATCC
 AGCTGCCTCAGATTGTCGTAGTGGAACTCAGAGCAGTGGAAAGAGCTCAGTGCTGGAAAGCCTAGTGGG
 GAGGGACCTTCTCCAGAGGGACTGGTGTGGTACCCGAGACCTCTTATTCTGCAGCTGGTCCACGTT
 TCACCAGAAGATAAAAGAAAAACAACAGGAGAAGAAAATGGAGTTGAAGCAGAAGAATGGGGTAAATTTT
 TTCACACCAAAAAACAAGCTTTACACAGATTTTGATGAAATTCGACAAGAAATGAAAATGAAACTGAAAG
 AATTTTCAGGAAATAATAAGGGAGTAAGCCCTGAGCCAATCCATCTCAAGGTTTTCTCGCCCAACGTTGTC
 AACCTGACGCTTGTGGATTTACCAGGAATGACCAAAGTACCTGTAGGTGATCAGCCCAAGGATATTGAGC
 TTCAAATCAGAGAACTTATTCTCGTTTATCAGTAATCCCAATCCATTATCCTTGCTGCTACTGCTGC
 AAATACAGATATGGCAACGTCAGAGGCACTTAAGATTTCCAGAGAGGTAGATCCAGATGGGCGCAGA
 ACTCTAGCTGTAATACCAAGCTTGATCTCATGGATGCGGGTACTGATGCCATGGATGTATTGATGGGAAGGG
 TTATTCCAGTCAAGCTTGAATAATTGGAGTAGTTAACAGGAGCCAACCTGGACATTAACAATAAGAAGAG
 TGTAAGTATTCAATCCGTGATGAGTATGCTTTTCTTCAAAGAAGTACCCATCTCTGGCCAACAGAAAT
 GGAACAAAGTATCTTGCTAGGACTCTGAATAGGTTACTTATGCATCATATCAGAGACTGTTTACCAGAGC
 TGAAAACAAGAATAAATGTCTTAGCTGCTCAGTATCAGTCTCTTCAAATAGCTATGGTGAGCCCGTGG
 TGATAAAAGTGTACTTTACTCCAGCTCATTACCAAGTTGCCACAGAGTATTGTAACACTATTGAAGGA
 ACGGCAAAGTACATTGAACTTCTGAGCTATGTGGTGGTGTAGGATTTGTTATATTTCCATGAGACTT
 TTGGGCGAACCTTAGAATCTGTTGATCCACTAGGTGGCCTAACACTATTGACATCTTGACCCGCAATTAG
 AAATGCTACTGGCCCCGTCAGCTTTATTTGTGCCTGAAGTTTCATTTGAGTTACTGGTCAAACGTCAG
 ATTAAGCGTCTAGAGGAGCCAGCTCCGCTGTGTGGAGCTGGTCCACGAGGAAATGCAGAGGATCATT
 AGCACTGTAGCAACTACAGCACACAGGAAGTACGGTTCCTAACTCCATGATGCCATAGTTGAAGT
 AGTGACTTGTCTTCTTCTGAAAAGGTTGCCCGTGACAAATGAAATGGTACATAACTTAGTGGCAATTGAG
 CTAGCATATCAACACAAAACACCCCGACTTTGCTGATGCTGTGGGCTAATGAACAATAACATAGAGG
 AACAAAGAAGAACAGACTAGCAAGAGAACTACCTTCCGCTGTATCGCGAGACAAGTTAATTCAGGACAA
 CAGAAGAGAACTAAAAATGTTGCATCTGCAGGTGGTGGAAATGGAGATGGTGGTGAATTGGAGATGGT
 GGTGAGGAAACGACAACAGGCAACTGGAGAGGAATGCTGAAAACCTCAAAGCTGAAGAATTACTTGCTG
 AAGAAAAATCAAACCAATTCCAATTATGCCAGCAAGTCCACAAAAGGCCATGCTGTCAATTTGCTAGA
 TGTGCCAGTTCAGTTGCACGAAAACCTGTGTCGCGGAGAACAGCGAGATTGTGAGGTTATTGAAAGACTC
 ATCAAATCATATTTTCTAATTGTGAGAAAGAAATATTCAAGACAGCGTCCCAAAGGCGAGTAATGCATTTT
 TGGTTAATCATGTGAAAGATACTCTCAGAGTGAATTAGTAGGGCAGCTGTATAAGTCACTTATTAGA
 TGACCTTCTGACTGAATCCGAGGACATGGCACAGCGAAGAAAAGAAGCGGCTGACATGCTGAAGGCATTA
 CAAGGAGCCAGTCAAATATTGCTGAAATCCGAGAGACTCATCTTTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_053655.3 , NP_446107.2
RefSeq Size:	2515 bp
RefSeq ORF:	2151 bp
Locus ID:	114114
UniProt ID:	O35303
Cytogenetics:	11q23
MW:	80 kDa
Gene Summary:	human homolog is a GTPase that is necessary for peroxisome division [RGD, Feb 2006]

Product images:



Circular map for RR205234