

Product datasheet for **RR205234L3V**

Dnm1l (NM_053655) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dnm1l (NM_053655) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dnm1l
Synonyms:	DLP1; Dnml1; Drp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053655
ORF Size:	2148 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205234).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053655.3 , NP_446107.2
RefSeq Size:	2515 bp
RefSeq ORF:	2151 bp
Locus ID:	114114
UniProt ID:	O35303
Cytogenetics:	11q23
Gene Summary:	human homolog is a GTPase that is necessary for peroxisome division [RGD, Feb 2006]



[View online »](#)