

Product datasheet for **RR205249**

Slc8b1 (NM_001017488) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc8b1 (NM_001017488) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc8b1
Synonyms:	Slc24a6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR205249 representing NM_001017488
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCAGGCAGATGGCTGGATCCGCTCTGGGCGCCTGGCTTCCTCTGTGTGGCTCTGATACTGGAGACGG
 CATCTGGGGCCGAGATCTATCCACAAAAGCACATGGACACATCCAGTTTTTCAGCCAGAGGTGTCAACCA
 GACTGCCATGGCAGACTGCCGAGCTGTGTGTAGCCTGAACACATCTGACCCTGTGACTTTGTCAAGAGG
 AATCCCGACTGCCACAGCGAGGGGGGCTACCTGGACTACCTCAAGGGCATCTTCTGCTACTTCCCCC
 ACCTCCTCCCTCTGGCCATCACCTCTATGTTTTCTGGCTCCTTTACCTCTTTCTGATCCTCGGAGTCA
 CGCAGCCAAGTCTTCTGCCCTAATTGTCAGCCATCTCCACCAGCCTCAAGCTCTCCACAACGTGGCA
 GGTGTACCTTCTGGCCTTTGAAATGGCGCTCCAGACATTTTCAGTGTCTCGTGGCTTTCTCAGACC
 CTCGACTGCTGGCCTGGCCATTGGGGCTCTGTTTGGTGCAGGGGTGCTGGTCACTACTGTGGTGGCTGG
 TGGCATACCATCCTGCGTCCCTCATGGCTGCCTCCAGGCCCTCCTCAGGGACATCACTTTCTATATG
 GTGGCCGTGTTCCCTCACCTTACCAGCCTCTATCTCGGCAGGATCACACTGGTGTGGGCACTGGGTTACC
 TGGGCCTCTATGTCTTCTATGTGGTACCGTCATCATCTGCACCTGGGTCTACCAAGGCAGCGGAGCAG
 GTCTTTGGTCCACTCCATATCAGAGACACCAGAGCTACTGACTGATTAGAAAGAGGATCAGATGTCTTCC
 AACACCAACAGCTATGACTATGGAGAGGAGTACCGGCCTCTGCTGCTGGGTGAGGAGACCACTGGGCA
 TCCTGCTCCAAGCCCTGAATCCCTTGGACTACAGGAAGTGGAGAAGTCACTCGATATCCTGCAAGCTCCT
 GAAAGTGGCCAAGTTGCCTGTGGAGTCTTGTGTGCTCACGGTACCGGTTGTGGACCCCGACAAGGAT
 GATCGGAACTGGAACGGCCACTCACTGTCTGCAGCTGGTCAACAGCCCTGGTCTGGTCTGACCC
 TGCAGTCGGGGTCTATGGCATCTATGAGATCGTGGCCTCCTTCTGTGTGGGCTGTAGTGGTATCGT
 GGGCACAGCCCTGGCTTCAAGTACCTTCTTGGCCACATCTAACAGCGAACCCCTAGACTGCACTGGCTC
 TTTGCTTCTGGGTTTCTGACCAGCGCCTGTGGATCAATGCAGCTGCCACGGAGGTGGTGAACATCT
 TGCGGTCCCTGGGCGTGGTCTTCCGCTCTGAGCAACACCGTCTAGGGCTGACCCTCCTGGCCTGGGGAA
 CAGCATTGGAGATGCGTTCTCAGATTTACGCTGGCCCGCCAGGGATACCCTCGGATGGCCTTCTCTGCC
 TGTTTTGGCGGCATCATTTCAATATCCTGTTGGTGTGGGACTGGGCTGCTGTACAAATGTCCGGA
 GCCATGCTCCGAGGTGAAGCTGGAACAGACGGACTGCTGGTGTGGGTGCTGGCCAGTGCCTGGGCT
 CAGCCTGGTCTTCTCCCTGGTCTCGGTTCCGCTCCAGTGTTCAGCTCAGCAAGGCCTACGGCCTCTGC
 CTCTCTCTTCTATATATGTTTATTGGTGGTCTCCTCACGGAGTTCGGGGTGATTCACTTGAAG
 CTGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR205249 representing NM_001017488
 Red=Cloning site Green=Tags(s)


MAGRWLDPLWAPGFLCVALIETASGAGDLSTKAHGHIQFSARGVNTAMADCRAVCSLNTSDRCDFVKR
 NPDCHSEGGYLDYKGFYFPPNLLPLAITLYVFWLLYLFLILGVTAAKFFCPNLSAISTSLKLSHNVA
 GVTFLAFNGAPDIFSAVAFSDPRTAGLAIGALFGAGVLVTTVVAGGITILRPFMAASRPFLRDITFYM
 VAVFLTFTALYLGRITLVWALGYLGLYVYVVTVICTWVYQRQSRSLVHSISETPELLTDSEEDQMSS
 NTSYDYGEEYRPLLLGEETTQILLQALNPLDYRKWRTQSI SCKLLKVAKLPVEFLLLL TVPVVDPDKD
 DRNWKRPLNCLQLVISPLVVLVTLQSGVYGIYEIGLLPVWAVVIVGTALASVTFATSNSEPPRLHWL
 FAFLGFLTALWINAATEVVNILRSLGVVFRLSNTVLGLTLLAWGNSIGDAFSDFTLARQGYPRMAFSA
 CFGGIIIFNILVGVLGCLLQIVRSHASEVKLEPDGLLVWVWLASALGLSLVFSVSVPLQCFQLSKAYGLC
 LLLFYICFIVVLLTEFGVIHLKAD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001017488

ORF Size: 1755 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017488.1](#), [NP_001017488.1](#)

RefSeq Size: 2693 bp

RefSeq ORF: 1758 bp


Locus ID: 498185

UniProt ID: [Q6AXS0](#)

Cytogenetics: 12q16
MW: 64 kDa
Gene Summary:

Mitochondrial sodium/calcium antiporter that mediates sodium-dependent calcium efflux from mitochondrion, by mediating the exchange of 3 sodium ions per 1 calcium ion. Plays a central role in mitochondrial calcium homeostasis by mediating mitochondrial calcium extrusion: calcium efflux is essential for mitochondrial function and cell survival, notably in cardiomyocytes (By similarity). Regulates rates of glucose-dependent insulin secretion in pancreatic beta-cells during the first phase of insulin secretion: acts by mediating efflux of calcium from mitochondrion, thereby affecting cytoplasmic calcium responses. Required for store-operated Ca(2+) entry (SOCE) and Ca(2+) release-activated Ca(2+) (CRAC) channel regulation: sodium transport by SLC8B1 leads to promote calcium-shuttling that modulates mitochondrial redox status, thereby regulating SOCE activity (By similarity). Involved in B-lymphocyte chemotaxis (By similarity). Able to transport Ca(2+) in exchange of either Li(+) or Na(+), explaining how Li(+) catalyzes Ca(2+) exchange. In contrast to other members of the family its function is independent of K(+) (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RR205249