

## Product datasheet for **RR205316L2V**

### Khdrbs3 (NM\_022249) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Khdrbs3 (NM_022249) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Khdrbs3
Synonyms:	Etle; SIm2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_022249
ORF Size:	1038 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205316).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022249.1</a> , <a href="#">NP_071585.1</a>
RefSeq Size:	1704 bp
RefSeq ORF:	1041 bp
Locus ID:	64015
UniProt ID:	<a href="#">Q9JLP1</a>
Cytogenetics:	7q34
Gene Summary:	mediates the selection of alternative splice sites; may provide a link between signal transduction pathways and pre-mRNA processing [RGD, Feb 2006]



[View online »](#)