

Product datasheet for **RR205337L4V**

Tmem176b (NM_134390) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmem176b (NM_134390) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tmem176b
Synonyms:	Lr8; Torid
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134390
ORF Size:	789 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205337).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134390.3 , NP_599217.1
RefSeq Size:	1683 bp
RefSeq ORF:	792 bp
Locus ID:	171411
UniProt ID:	Q925D4
Cytogenetics:	4q24
Gene Summary:	Required for the development of cerebellar granule cells (By similarity). May play a role in the process of maturation of dendritic cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)