

Product datasheet for **RR205348L4V**

Acrbp (NM_001025049) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acrbp (NM_001025049) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Acrbp
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025049
ORF Size:	948 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205348).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025049.1 , NP_001020220.1
RefSeq Size:	1284 bp
RefSeq ORF:	951 bp
Locus ID:	500316
UniProt ID:	Q6AY33
Cytogenetics:	4q42
Gene Summary:	May be involved in packaging and condensation of the acrosin zymogen in the acrosomal matrix via its association with proacrosin.[UniProtKB/Swiss-Prot Function]



[View online »](#)