

Product datasheet for **RR205355**

Nrn1l (NM_001135006) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nrn1l (NM_001135006) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nrn1l
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205355 representing NM_001135006
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGATGTGCTACTGCTGCCACTGCCACTGGCGCCCGCTGTGAGCGGTACCCTGTGCCCTGACTGTG
TGCTGCTACTGCCACTCGCAGTGGCCTCTGAGGGCCAAACCGCTGTGAGACCATATACCAAGGCTTTGC
TGAATGTCTCATCCGCTGGGGGATGGCATGGGCCGAGGAGGCGAGCTACAGACCGTCTGCAGATCCTGG
AATGAATTCACGCCTGTGCCTCTCGGGTCTGTGAGGCTGCCAGAGGAGGCTGCTGCAGTGTGGGAGT
CACTGCAGCAAGAAGCTCGCCGTGCCCCACACCCAGATAAATTGCACATCCTCTGTGGCGCTCCTGTGAG
TGTTCCGGAGCGGATTGCTGGCCAGAGACCAACCAGGAGACTACGGGCCACAGCTCCTGCCCTGGCT
TCAGCCCCAGCCCTGCGTTGCTCGCTGCTGCTGGCGCTTGCCTGCCTCCTGGGGCTCTGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR205355 representing NM_001135006
Red=Cloning site **Green**=Tags(s)

MMCYCCHCHWRRRCQRLPCALLLLLLPLAVASEGPNRCETIYQGFAECLIRLDGMGRGGELQTVCRSW
NEFHACASRVLSGCPPEAAAVWESLQEQEARRAPHDPNLHILCGAPVSVRERIAGPETNQETLRATAPALA
SAPAPALLAAALALACLGLPLA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001135006

ORF Size: 486 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001135006.1, NP_001128478.1</u>
RefSeq Size:	620 bp
RefSeq ORF:	489 bp
Locus ID:	100189544
Cytogenetics:	19q12
MW:	17.4 kDa