

Product datasheet for **RR205375**

Fzd2 (NM_172035) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fzd2 (NM_172035) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fzd2
Synonyms:	Fz-10; Fzd10
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR205375 representing NM_172035
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCGGGCCCGCAGCGCCCTGCCCGCAGCGCCCTGCCCGCCTGCTGCTGCCACTGCTGCTGCTGCCGG
 CTGCCGGCCCGCTCAGTTCCACGGGGAGAAGGCATCTCCATCCCGGACCACGGCTTCTGCCAGCCAT
 CTCCATCCCGCTGTGCACGGACATCGCCTACAACCAGACCATCATGCCAACCTTCTTGGGCACACGAAC
 CAAGAGGACGCGGGCCTGGAGGTGCATCAATTCTACCCGCTGGTGAAGGTGCAGTCTCGCCGAGCTGC
 GCTTCTTCTGTGCTCCATGTACGCTCCGGTGTGCACGGTGTGGAGCAGGCCATCCCGCGTCCGCTC
 CATCTGCAACCGCGCGCCAAGGCTGCGAGGCGCTCATGAACAAGTTCGGTTCAGTGGCCCGAGCGC
 CTCGCTGCGAGCATTTCCCGCTCACGGCGCGGAGCAGATCTGCGTGGCCAGAACCACTCCGAGGACG
 GAACTCCTGCGCTACTACCACCGCGCCACCGTCTGGGCTGCAGCCTGGCGTGGTGGCACCCCGGGCGG
 CCTGGCGGTGGTGGCGGCCCGCGCTACGCCACTCTGGAGCACCTTTCAGTGTCCCGCGTCTCT
 AAGGTGCCGTCTATCTCAGCTATAAGTTTCTGGGTGAGCGGATTTGTGCCGCGCCCTGCGAGCCTGCAC
 GGCCCGACGGCTCCATGTTCTTCTCGCACCCACCACTCGTTTTGCCCGTCTCTGGATCCTCACATGGTC
 GGTGCTGTGCTGCGCTTCTACTTTCTTACGGTACCACCTCTTAGTGGCCATGCAGCGATTCCGCTAC
 CCAGAGCGGCCATCATCTTCTGTCCGGTGTACACCATGGTGTGAGTGGCCATACATTGCGGGCTTCG
 TGCTCCAGGAGCGGTGGTGTGAACGAGCGTCTCTGAGGACGGTTATCGCACGGTGGGGCAGGGCAC
 TAAGAAAGAAGGCTGTACTATACTTTCATGATGCTCTACTTCTCAGTATGGCCAGTCCATCTGGTGG
 GTGATTCTGTCCCTCACCTGGTCTCTGGCAGCCGGTATGAAGTGGGGCCACGCGCCATCGAGGCCAATT
 CGCAGTACTTCCACCTGGCCGCTGGCGGTGCCCGCCGTCAAAACCATCACCATCTCGCCATGGGCCA
 GATCGACGGCGACCTGCTGAGCGCGTGTGCTTCTGTTGGCCTCAACAGGCTGGACCCGCTGCGAGGCTTC
 GTGCTGGCGCCGCTCTTCTGTACCTGTTTCATCGGCACATCCTTCTGCTGGCGGGCTTCTGTACTCT
 TCCGATCCCGACCATCATGAAGCACGACGGCACCAAGACGGAGCCGCTGGAGAGGCTCATGGTGCAT
 CGGCGTCTTCTCCGTGCTTACACCGTACCAGCCACCATCGTATCGCCTGCTACTTCTATGAGCAGGCC
 TTCCGCGAGCACTGGGAGCGCTCGTGGTAAGCCAGCACTGCAAGAGCCTAGCCATCCCTGCCCGGCC
 ACTACACGCCCGCACGTCGCCGACTTACAGTCTACATGATCAAATACCTCATGACGCTCATCGTGGG
 CATCACGTCCGGCTTCTGGATCTGGTCCGGCAAGACGCTGCACCTCGTGGAGGAAGTCTACACCGCTCTC
 ACCAACAGCCGGCATGGAGAGACCACCGTG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTAA

Protein Sequence:

>RR205375 representing NM_172035
 Red=Cloning site Green=Tags(s)


MRARSALPRALSALPRLLLPLLLLPAAGPAQFHGEKGISIPDHGFCQPISIPLCTDIAYNQTIMPNLLGHTN
 QEDAGLEVHQFYPLVKVQCSPFLRFLCSMYAPVCTVLEQAIPPCRISICERARQCEALMNKFGQWPER
 LRCEHFPRHGAEQICVGNHSEDTGTPALLTTAPPSTGLQPGAGGTPGGPGGGGAPPRYATLEHPFHCPVRL
 KVPSYLSYKFLGERDCAAPCEPARPDGSMFFSHHTRFARLWILTWSVLCASFFVTVTTSLVAMQRFYR
 PERPIIFLSGCTMVSVAIYAGFVLQERVVCNERFSEDGYRTVGQGTKEGCTILFMMLYFFSMASSIWW
 VILSLTWFLAAGMKWGHAAIEANSQYFHAAWAVPAVKITITILAMQIDGDLISGVCVGLNRLDPLRGF
 VLAPLFFVYLFIGTSFLLAGFVSLFRIRTIMKHDGKTEPLERLMVRIGVFSVLYTVPATIVIACYFYEQA
 FREHWERSWVSQHCKSLAIPCAHYTPRTSPDFTVYMIKYLMTLIVGITSGFWIWSGKTLHSWRKFYTRL
 TNSRHGETTV

SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_172035
 ORF Size: 1710 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172035.1 , NP_742032.1
RefSeq Size:	1912 bp
RefSeq ORF:	1713 bp
Locus ID:	64512
UniProt ID:	Q08464
Cytogenetics:	10q32.1
MW:	63.9 kDa
Gene Summary:	G-protein coupled receptor that binds Wnt proteins and induces intracellular calcium release; also mediates signaling via inhibition of cGMP [RGD, Feb 2006]