

## Product datasheet for **RR205386L3V**

### Smagp (NM\_182817) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Smagp (NM_182817) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Smagp
Synonyms:	MGC124543
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_182817
ORF Size:	291 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205386).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182817.1</a> , <a href="#">NP_877969.1</a>
RefSeq Size:	1501 bp
RefSeq ORF:	294 bp
Locus ID:	300236
UniProt ID:	<a href="#">Q7TPF1</a>
Cytogenetics:	7q36
Gene Summary:	May play a role in epithelial cell-cell contacts. May play a role in tumor invasiveness and metastasis formation.[UniProtKB/Swiss-Prot Function]



[View online »](#)