

Product datasheet for **RR205436**

Il17re (NM_001004091) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Il17re (NM_001004091) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Il17re
Synonyms:	MGC93201
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR205436 representing NM_001004091
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGAGCCCCAGACTGGCAGCCTTGCTCCTGTCTCAACCTCTACTGTTTCATCTGCCTTGCTGTCTG
CTCAGGTTGCCTGCCCTGCCTGCCACGTTGGACCAGCCACTGTCTACTGGCCTTCGGTGTGGATAAGCA
TTTTGCTGGCCTTCAGTGTGGCTGGTCCCTCTCTTGGTGAGGAAATCTAAAAGTTCTCATGAATTTGAA
GTCTGTTGGAGGCACTGGACACCATCACCTTCCAGAGGAAGCTACTAGGTAGCCCTTCCTGTCTCAGG
AGAGCCATCGAATTTCCACCCACTTCTAGCCATCTCCACAGAGGCTACGTACCAAAAGGACCCAGCC
TTCAGCTGCAGAAGGGATAGAACATCTCCCTGGAGCAGGTTACAAAAAGCATGGAGGACCTGAATTTCC
TTTGATTTGCTGCCTGAGGTACAGGCTATTCGGGTGACTATTCCTCCAGGCCAGAGGCCAACGTACGCC
TTTGTACCAGTGGGCACTGGAATGTGAGGACATGAGTAGTCCTTTTGATACCCAGAAAATTGTGTCTGG
GGCCCACTGTAGACCTGCCTTATGAATTCCTTCTGCCTTGATGTGCATAGAGGCTCCTACCTGCAG
GAGGACACTGTGAGGCACAAAAAGTGTCCCTTCGGAGCTGGCCTGAAGCTTATGGCTCCGACTTCTGGA
AGTCAATACGCTTCACTGACTACAGCCAGCACAGTCAGATGGTCATGGCTCTGACTCCGCTGCCCACT
GAAACTGGAGGCCTCTCTGCTGGAGGCAGGCCAGTACGCCCTGCAAAATCCTTCCCAATGCCACA
GCAAAGGAGTCAGAAGGATGGTATATCCTGGAGAATGTGGACTTGACCCCCAGCTCTGCTTTAAGTTTT
CTTTTGAAAACAGTAGCCATGTTGAATGTCCCCACAGAGTGGCTCTCTCCATCCTGGACTGTGAGCAT
GGATACCCAGGTTTGGGCTGGATAGCTCCATGACTCCTGTGTACAGCATCAGTCAGACCCAGGGCTCAG
TTCCAGTGACACTGAACCTCATTATTCCTTCTTGAGGCAGGGGAGCTGCATCCTGGTGTGGAGGTCAGA
TGCCAGTTTGCCTGGAAGCACCTCTTGTGTCTGATGTCTCCATAGACACCTCGGGCTCTTGATCCTG
GCACTGCTGGGTCTACCACTCTACTGGGTGATGTTCTGGTCTCTTCTGCCGGCCTACTGCCAGGCC
CCGGTCGAACAAGGCCAGTTTTACTCTACACGCAGCGGACTCAGAGGCACAGCGACGCTGGTGGGAGC
TTTGGCCGAAGTCTGCGGACCGCCTGGGTGGTGGACGCGACGTGATCGTGGATCTCTGGGAAGGGACG
CACGTGGCTCGCATTGGACCACTGCCATGGCTCTGGGCAGCACGGGAACGCGTGGCACGGGAGCAGGGCA
CTGTGTTGCTCCTTTGGAGCTGTGCGGGCCCCAGCACTGCCTGCAGCGGCGACCCGAGACTGCATCCCT
TCGCACCTTGTCGTGCGCTGCTCCACGCCAGCTGCTGCTCGCCTACTTCAGTCGCCTCTGTGCCAAAGT
GACATTCGGGGCCGCTGCGCGCTACCGCGCTACCGCCTGCTTCGTGATCTGCCACGCCTACTGAGAG
CGCTGGATGCTCGACCTGCCACCCTAGCCACTAGCTGGAGTACCTTGGAGCTAAGGAGTGTGAAAAG
TCGCCTGGAGCTGTGTCACCTGCTGGAACCTCGAGGCTGCCAAAGACGACTACCACGGTTCACCAATAGT
CCTTGTGGCTTCAGTTGTCTG

AC**CGGGCCGC**TCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA
TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR205436 representing NM_001004091
 Red=Cloning site Green=Tags(s)


MGSPRLAALLLSQPLLFICLAVSAQVACPCLPRWTSCHLLAFRVDKHFAGLQCGWFPLLVRKSKSSHEFE
 VCWRHWTPSPFQRKLLGSPSLSQESHRISTHFLAISHRGLRTRKRTQPSAAEGIEHLPGAGSQKHGGPEFS
 FDLLPEVQAIRVTIPPGPEANVRLCYQWALECEDMSSPFDTKIIVSGGHTVDLPYEFLLPCMCI EASYLQ
 EDTVRHKKCFRSPWPEAYGSDFWKSIRFTDYSQHSQMVMLTLRCPLKLEASLCWRQDPVTPCKILPNAT
 AKESGWYILENVDLHPQLCFKFSFENSSHVECPHQSGSLPSWTVSMDTQAQQLILHFSTRTYATFSAAW
 SNPGLGLDSSMTPVYSISQTQGSVPVTLNLIIPFLRQGCILVWRSDVQFAWKHLLCPDVSHRHLGLLIL
 ALLGLTLLGVVLVLCRRLLPGRTRPVLLHAADSEAQRRLVGALAE LLRTALGGGRDIVDLWEGT
 HVARIGPLPWLWAARERVAREQGTVLLWSCAGPSTACSGDPQTASLRTLSCAAPRQLLLAYFSRLCAKG
 DIPGPLRALPRYRLRLDLPRLLRALDARPATLATSWSHLGAKECLKSRLELCHLLEEA AKDDYHGSTNS
 PCGFSC L

TRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-NotI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:


ACCN:	NM_001004091
ORF Size:	1911 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001004091.1 , NP_001004091.1

RefSeq Size: 2640 bp

RefSeq ORF: 1914 bp

Locus ID: 362417

UniProt ID: [Q6AZ51](#)

Cytogenetics: 4q42

MW: 70.9 kDa

Gene Summary: Specific functional receptor for IL17C, signaling through the NF-kappa-B and MAPK pathways. Requires TRAF3IP2 /ACT1 for signaling. Crucial regulator in innate immunity to bacterial pathogens (By similarity).[UniProtKB/Swiss-Prot Function]