

## Product datasheet for RR205445L3

# Rasgrp1 (NM\_019211) Rat Tagged Lenti ORF Clone

### **Product data:**

**Product Type: Expression Plasmids** 

**Product Name:** Rasgrp1 (NM\_019211) Rat Tagged Lenti ORF Clone

Tag: Myc-DDK Symbol: Rasgrp1 Synonyms: Rasgrp **Mammalian Cell** 

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Puromycin

E. coli Selection: Chloramphenicol (34 ug/mL)


The ORF insert of this clone is exactly the same as(RR205445). **ORF Nucleotide** 

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn



#### Plasmid Map:



**ACCN:** NM\_019211 **ORF Size:** 2385 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 019211.2</u>, <u>NP 062084.2</u>

RefSeq Size: 3624 bp RefSeq ORF: 2388 bp



### Rasgrp1 (NM\_019211) Rat Tagged Lenti ORF Clone - RR205445L3

 Locus ID:
 29434

 UniProt ID:
 Q9R1K8

 Cytogenetics:
 3q35

**Gene Summary:** putative guanine nucleotide exchange factor with Ca2+ and diacylglycerol binding domains;

may activate Rap protein mediated inhibition of Ras signaling pathways [RGD, Feb 2006]