

Product datasheet for **RR205484**

Bhlha15 (NM_012863) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bhlha15 (NM_012863) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bhlha15
Synonyms: bHIH; Bhlhb8; Mist1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205484 representing NM_012863
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGACCAAGAACCGGCCCTCGGCGCCGGACACCCATGCAGGACGCAGAAGCCACCCAGGAGAGC
 AGACGCCTGACAGATCCCAGTCAGGCTCAGGGGCATCAGAAGTGACCAAGGGTCTTCGGAGCAGGACAGC
 CCGTGCAGCGGAACTCGGGCAGAGGTCAGCCCGCGGACAGGGGTCCAGCAGCCGAGGAGAACAGC
 GTTCAGAGGGCGGCTGGAGAGCAATGAGCGAGAAAGACAACGGATGCACAACTGAACAATGCTTTCCAGG
 CGCTGCGGGAGGTCATCCCGCATGTGCGGGCTGACAAGAAGCTCTCCAAGATCGAGACCCCTCACACTGGC
 CAAGAATAATCAAGTCGTCGACAGCCACCATACTTACTATGTCCAGCAGCCGCCTCCCGGGGCTGGAG
 GCGCCAGGTCCTGCGCCAGGCCCTAAATTGTACCAGCATTACCATCACCAGCAGCAGCAACAGCAGCAGC
 AGCAGCAGGTGGCTGGGGCCGTGCTCGGTGTCAACGAGGACCAGCCTCAAGGCCACCTGCAACGATACTC
 TACACAGATCCACAGCTTCAGAGAGGGGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR205484 representing NM_012863
 Red=Cloning site Green=Tags(s)

MKTKNRP RRRTPMQDAEATPGEQTPDRSQSGSGASEVTKGLRSRTARASGTAEVSRRRQGSRRRENS
 VQRRLESNERERQRMHKLNNAFQALREVIPHVRADKLSKIETLTLAKNYIKSLTATILTMSSSRLPGLE
 APGPAPGPKLYQHYYHQQQQQQQVAGAVLGVTEDPQPQGLQRYSTQIHSFREGS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_012863

ORF Size: 591 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_012863.1](#), [NP_036995.1](#)

RefSeq Size: 3546 bp

RefSeq ORF: 594 bp

Locus ID: 25334

UniProt ID: [P70562](#)

Cytogenetics: 12p11

MW: 22.2 kDa

Gene Summary: basic helix-loop-helix transcription factor; may regulate gene expression in embryonic and postnatal cell lineages [RGD, Feb 2006]