

## Product datasheet for **RR205494**

### Ankrd6 (NM\_001134969) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ankrd6 (NM_001134969) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ankrd6
Synonyms:	MGC188478
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RR205494 representing NM\_001134969  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGCCAGCAGGATGCGGTGCGCCGCTTTTCAGAACGCCTGCTCATAGCTGCGTACAAAGGCCAAACAG  
AGAATGTGGTTTCAGTCACTCAACAAGGGTGCCAAGGTAGCGGTTACCAAGCATGGCCGGACCCCTTGCA  
CCTCGCTGCCAATAAGGGTCATCTTTCTGTGGTGCAGATTTTGCTGAAGGCTGGCTGTGACCTCGATGTA  
CAAGATGATGGGACCAGACAGCCTTGACCCGGCCACTGTGGTGGGCAACACTGAGATCCTCACAGCGC  
TTATCCGGGAGGGATGTGCCCTGGACAGACAGGACAAGGATGGGAATACGGCCCTGCATGAAGCCGCTG  
GCATGGATTACGCCAGTCAGCCAAGCTGCTCGTTAAAGCCGGAGCCAACGTGCTTGCCAGAAAACAAGGCA  
GGGAACACGGCTCTGCACCTGGCCTGCCAGAATAGTCACTCCCAGAGCACACGCGTCTCTCTGTTGGCG  
GCTCCCGGGCTGACCTCAAAAAAATGCAGGGGACACCTGCTTGACGTGGCTGCGCGCTATAACCACTT  
GTCCGTTGTTAGACTCCTCCTCAGTGTCTTCTGCTCTGTTTCATGAAAAGAACCAGGCTGGAGACACGGCC  
CTCCATGTGGCCGCTGCCCTAAATCACAGAAGGTAGTTAAAAATCCTGCTGGAAGCCGGTCTGACACGA  
CCCTTGTGAACAATGCAGGGCAGACCCCACTGGAGACCCCGCTACCACGCAATCCCAGGATTGCTCT  
GCTCCTACCAAGGCTCCCAGATCTTGCCTTCAGTCGTGGCGAAGCCTGAGGAAAAGGAGGGAGAGG  
CTCAAGGAAGAGAGGAGGGCTCAGTCTGTGCCAGAGATGAGGTGGCCAGAGCAAGGGAAGCGTCTCAG  
CAGGAGACACCCCAAGCAGTGAACAAGCTGTACCCGAAAAGGAAGAAGCTAGGAGAGATTGCCACCAGC  
TTCCCAGGAGCCAAGAAAGGATGAAAGGAGGAGAAAGTCAAGGCCAGAGGTGTCAGCACTATCTGACCC  
ACCCAGCAGCAGACCAACAACCTGGACACCAAAAGAACCTACATGGTCATCACCACCAAAAAAGAAGA  
GCAGACATCGATGCTGGTCCCAACCCCGGGTTCAGGGCGTACCAGCTCTACACGCTGTACCGTACCG  
TGGCGAGGACGGAAAGTGATGCAGGCACCAATCAAAGGTTGTCGCTGTGAGCCTCTGATCAACAAGCTG  
GAGAACCAGCTGGAGGCGCTGTGGAGGAGTTAAGAGCTGAGCTGGGAACAGTGCAGGACAAAGTGAACA  
CAAAGCTGGGCCACATGGAGAGCAAGAGCCAGCACCAATGTGCGTTTTGGACAACTGATGGTGGAGCG  
GCTGTCAGCTGAGAGAACGGAATGTCTGAATCGCTGCAGCAGCATGCAGCTGCTGAGAAGCAGGAGGGA  
GAGAAACGTGAGATGTCATTGGTGGATGAACTGAAATCCTGGTGTATGTTGAAGATCCAGAATCTGGAGT  
TGAGACTTTCGGGAGAGTCTCGGACCTTAGAGCTAAATCCACACCATCTCCCTCCGACTCCTCTCCAGC  
TGTAGACCAACCAGTGGTAGCTGCAGGGCCGGGAGCAGCTCCGACAGCTCCTCTCAGGTTGTGAGGCC  
AAGGACAAGGCACTCAATGCCACTGCTGTCCCAAGACACCAGCAGGAGCTGCTTCCCTCTGATTGCACAG  
GCTCCGGGTTAAGGAAGGTCAAGGCCCGGCTGCTTCCAGACTGGGTGATCAGCAGACTGGGTCTGTGT  
CAACCGAGGCACTCAAACCAAGAAGTCTGCGAGAGGTGGGACAGCAAGCATCGAGGCCAACAACCGACA  
GCCAGCAGCCCAAGTGGGACAGCCATCCGCGCAAGCAGTGCAGCCCGAGATGCCTCGCAAGCCCTGG  
AGCTAACCCAGTATTTTTTTGAGGCTGTTTCTGCTCAGATGGAAAAGTGGTACGAACGGAAAATTGAAGA  
AGCAGCAAGCCAGGCCAGTCAGAAAGCTCAGCAAGATGAGGCCCACTGAAGGAGCACATCAGGAGTCTA  
GAAGAGGAGGTTGCCAAGCTGCGCACAAAGATGCACAAG

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA


**Protein Sequence:** >RR205494 representing NM\_001134969  
 Red=Cloning site Green=Tags(s)

```
MSQQDAVAALSERLLIAAYKQQTENVVQLINKGAKVAVTKHGRTPLHLAANKGHL SVVQILLKAGCDLDV
QDDGDQTLHRATVVGNTEIL TALIREGCALDRQDKGNTALHEAAWHGFSQSAKLLVKAGANVLARNKA
GNTALHLACQNSHSQSTRVLLLGGSRADLKNNAGDTC LHVAARYNHL SVVRLLLSAFCSVHEKNQAGDTA
LHVAAAALNHKKVVKILLEAGADTTLVNNAGQTPLETARYHDNPEVALLLTKAPQILRF SRGRSLRKRER
LKEERRAQSVPRDEVAQSKGSVSAGDTPSEQAVPEKEEARDCPPASQEPKDERRRKSRPEVSALSDP
TPAADQQPGHQKNLHGHHHQKKKSRHRCWSPPPHGFRA YQLYTLYRGEDGKVMQAPIKGCRCPELINKL
ENQLEAAVEELRAELGTVDKVN TKLGHMESKSQHQM CVL DKL MVERLSAERTECLNRLQQHAAAEKQEG
EKRQMSLVDELKSWMLKI QNLELRLSGESRTFRAKSTPSPSDSPA V DQPVVAAGPAAASDSSSQVVRP
KDKALNATAVPRHQE L LPSDCTGSLRKYKAPAASRLGDQQTGSCVNRGTQTKKSARGGQTKHRGQQPT
ASSPSGQQPSAASSDARDASQALEL TQYFFEAVSAQMEKWYERKIEEAR SQASQKAQQDEATLKEHIRSL
EEEVAKLR TKMHK
```

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_001134969

**ORF Size:** 2139 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001134969.1](#), [NP\\_001128441.1](#)

**RefSeq Size:** 4638 bp


**RefSeq ORF:** 2142 bp

**Locus ID:** 500430

**Cytogenetics:** 5q21

**MW:** 78.1 kDa

**Product images:**



Circular map for RR205494