

## Product datasheet for **RR205535L4V**

### Zbtb10 (NM\_024489) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zbtb10 (NM_024489) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zbtb10
Synonyms:	Rinzf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024489
ORF Size:	2508 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024489.1</a> , <a href="#">NP_077815.1</a>
RefSeq Size:	3849 bp
RefSeq ORF:	2511 bp
Locus ID:	80338
UniProt ID:	<a href="#">Q9WTY8</a>
Cytogenetics:	2q23
Gene Summary:	binds the gastrin CACC cis-regulatory element; may block Sp1 transactivation of gastrin gene expression [RGD, Feb 2006]



[View online »](#)