

Product datasheet for **RR205537**

Rnase11 (NM_001012476) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rnase11 (NM_001012476) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rnase11
Synonyms: Raj1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205537 representing NM_001012476
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCAGCCTTTCTGCTGCTGCTTGTCTTGGACTGCTTCTTGTAGAGCCTTCGAAAGCAAACGAAGG
 GGGCCAGGGAACAATTTCCAGGAAGAGACTCGATCTGCTGCAAACAGACTCCTGAGGAATCAACAAA
 CTC AACCTGT CAGATGAAAACATCAGCCTCGGCATATCCAGGCACGTGATGTCTGCCCCACCCGAACA
 TC CAGAAGGCTGTCTTTTGT CATCCCCAAGAGAAATACTGTGAGGAATGGCAGGGACTGTGTAATAGCC
 CGAGTGTGTGGAGAACAGATGTGGAGGTGAATGAGTCTTGCCAGTTGGCAATAACTTTATACATAGCTC
 CATGGATGTGAGCCACGGAATCCCAAGGCCACCAGCGGGAAGTGCGAACGAACACCTAACCTAAGATAC
 TGTGATAGCTTAGGACTGGAGTGTACCATGTGCAAGGTCCTTGCAAGGCCACCAGTGCCCCAGGTACCATG
 AGCACAGGATAACCTCATTAAAGAGAATCTTGACGGTGTGACAAGCCATTCTCTGATGAGCTGTTAGT
 TACTGGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR205537 representing NM_001012476
 Red=Cloning site Green=Tags(s)

MAAFLLLLVLG LLLVEPSESKTKGAREQFSQEETRSAAKQTPEESTNSTLSDENISLGISRHVMSAPPRT
 SRRLSFVIPKRNTVRNGRDCVNSPSVWR TDVEVNESQLGNNFIHSSMDVSHGIPKATSGKCERTPNLRY
 CDSLGLLECTMCKVLGHQCPRYHEHRITSLKRILTVLTSLSLMSWLVTGC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001012476

ORF Size: 570 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001012476.2](#), [NP_001012494.1](#)

RefSeq Size: 1503 bp

RefSeq ORF: 573 bp

Locus ID: 497651

Cytogenetics: 15p14

MW: 21.1 kDa