

Product datasheet for **RR205544L3V**

Stx2 (NM_012748) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stx2 (NM_012748) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Stx2
Synonyms:	Epim
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012748
ORF Size:	870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205544).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012748.1 , NP_036880.1
RefSeq Size:	911 bp
RefSeq ORF:	873 bp
Locus ID:	25130
UniProt ID:	P50279
Cytogenetics:	12q13
Gene Summary:	Essential for epithelial morphogenesis. May mediate Ca(2+)-regulation of exocytosis acrosomal reaction in sperm.[UniProtKB/Swiss-Prot Function]



[View online »](#)