

Product datasheet for **RR205547**

Casp2 (NM_022522) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Casp2 (NM_022522) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Casp2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205547 representing NM_022522
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCGTCGAGCGGGAGGTCTCAGTCTTCCCTGCACAGGAAGGGGCTGATGGCGGCTGACAGGAGGA
GCAGGATTTTGGCAGTGTGTGGAATGCATCCTGACCATCAGGAACTCTGAAAAAGAATCGAGTGGTGCT
GGCCAAACAGCTGCTGCTGAGCGAACTGTTAGAACATCTCCTAGAAAAGGACATTATCAGTTTGGAAATG
AGGGAGCTAATCCAGGCCAAAGGGGCGAGTTTCAGCCAGAATGTGGAACCTCTCAACCTGCTGCCAAAGA
GAGGACCCAGGCTTTTGTGCCTTCTGTGAAGCCCTGCGGGAGACCAGGCAGGGTCACTTGAAGACTT
ACTGCTCACAAACCCTCTCTGATATCCAGCACATACTCCCACCGTTGAGCTGTGACTATGACTCAAGTCTC
CCCTTCTCGGTGTGTGAGTCTGCCCTCTCACAAAGCAGTCCCGCCTATCCACAGATACTATGGAACACT
CCTTAGATAATGGTGATGGTCCCTGTCTTCAGGTGAAGCCATGCACTCCTGAGTTTTACCAGGCACA
CTACCAGCTGGCCTATAGGTTGCAGTCTCAGCCTCGTGGCTTGGCACTGGTGATGAGTAATGTGCACTTC
ACTGGAGAGAAAGACCTGGAATTTGCTCCGGAGGGGATGTGGACCACACTACTAGTCAACCCTTTCA
AGCTTTTGGGCTACAATGTCCATGTGCTGTATGACCAGACTGCACAGGAAATGCAAGAGAAACTTCAGAA
TTTTGCACAGTTACCTGCACACCGGTCACAGACTCCTGCATAGTAGCACTCCTCTCACATGGTGTGGAA
GGTGGCATCTACGGGGTGGATGGCAAATGCTTACAGTCCAAGAGGTTTTTCGACTTTTTGACAATGCTA
ACTGTCCAAGTCTACAGAACAACCAAAAATGTTCTTCAATCCAAGCATGTCGTGGAGATGAGACAGATAG
AGGGGTTGACCAGCAAGATGGAAGAACCATGCACAATCCCTGGATGTGAGGAGAGTGTGCTGGCAAG
GAGGAATTAATGAAGATGAGACTTCTACTCGCTCAGACATGATATGTGGCTATGCTTGCCTAAAGGCA
ATGCTGCCATGCGGAATACCAACGGGGTTCCTGGTACATTGAGGCCCTCACTCAGGTGTTCTCTGAGAG
AGCTTGTGACATGCACGTGGCCGACATGCTTGTAAAGTGAATGCCCTTATCAAGGAGCGTGAAGGCTAT
GCCCTGGTACAGAATTCACCGATGCAAGGAGATGTCTGAGTACTGTAGTACTGTGCCAGCAGCTCT
ACCTGTTCCAGGCTACCCACCCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR205547 representing NM_022522
 Red=Cloning site Green=Tags(s)


MAASSGRSQSSLHRKGLMAADRRSRILAVCGMHPDHQETLKKNRVVLAKQLLSELLEHLLLEKDIITLEM
 RELIQAKGGFSQNVLLNLLPKRGPQAFDAFCEALRETRQGHLEDLLLTTLSDIQHILPPLSCDYDSSL
 PFSVCESCPCHKQSRLSTDTMEHSLDNGDGPCLQVKPCTPEFYQAHYQLAYRLQSQPRGLALVMSNVHF
 TGEKDLFRSGGDVDHTTLVTLFKLLGYNVHVL YDQTAQEMQEKLQNF AQLPAHRVTDSCIVALLSHGVE
 GGIYGVGKLLQLQEVFRLFDNANCP SLQNKPKMFF IQACRGDE TDRGVDQDQDGKNHAQSPGCEESDAGK
 EELMKMRLPTRSDMICGYACLKGNAAAMRNTRKGSWYIEALTQVFSERACDMHVADMLVKVNALIKEREGY
 APGTEFHRCKEMSEYCSLTCQQLYLFPGYPPT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_022522
ORF Size:	1356 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022522.2 , NP_071967.2
RefSeq Size:	3365 bp
RefSeq ORF:	1359 bp
Locus ID:	64314
UniProt ID:	P55215
Cytogenetics:	4q24
MW:	50.7 kDa
Gene Summary:	may play a role in proteolysis and programmed cell death [RGD, Feb 2006]