

Product datasheet for **RR205590L4V**

Eef1g (NM_001004223) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eef1g (NM_001004223) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Eef1g
Synonyms:	MGC94929; MGC114210
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001004223
ORF Size:	1311 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205590).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004223.2 , NP_001004223.1
RefSeq Size:	1433 bp
RefSeq ORF:	1314 bp
Locus ID:	293725
UniProt ID:	Q68FR6
Cytogenetics:	1q43
Gene Summary:	Probably plays a role in anchoring the complex to other cellular components. [UniProtKB/Swiss-Prot Function]



[View online »](#)